

**ACCURACY**

Pressure lock, which keeps the tube stabilized.




**HYGIENE**

With opening that allows hygiene and access to the oral cavity avoiding damage or injury to the patient.

**SAFETY**

Fixing band with Velcro adjustment, which ensures the correct positioning of the device.

TRACHEOFIX,  
THE  
TRACHEOTRACHEAL PLUS TUBE  
HOLDER,  
PROVIDES AN  
EFFECTIVE WAY  
TO KEEP THE  
TUBE SAFE AND  
COMFORTABLE  
ON THE PATIENT.



**COMFORT**  
Medical-grade silicone bite, ensuring  
greater patient comfort.

**TRACHEOFIX PLUS**

- INCLUDED ITEMS**
- Tracheofix + Plus
  - Velcro clamping band

Latch Material	Polycarbonate
Body Material	Polycarbonate
Band Material	Polyester/Nylon
Cushion Material	Silicone

**INDICATION**

It is intended to maintain safe fixation in intubated patients under mechanical ventilation, reducing the risk of accidental extubation.

551   TRACHEOFIX + GMI ENDOTRACHEAL TUBE CLAMP PLUS ANVISA registration No. 80423540067			
Code	BODY LENGTH (cm)	SILICONE BITE LENGTH (cm)	Size
551.1.1	09	2.1	Paediatric
551.3.1	10	3.1	Adult