


SAFETY

Protects and unblocks the airway against tongue falling over the posterior pharynx.
Protects the endoscopic probe from possible bite trauma.
Prevents dental trauma to the patient.



ACCURACY

Two-way PVC, responsible for the administration of O2, and capnography reading.




OROPHAR, THE
GMI DOUBLE-
LUMEN
OROPHARYNGEAL
CANNULA,
MAINTAINS THE
AIRWAY OPEN
DURING
ENDOSCOPIC
PROCEDURES,
BRINGING SAFETY
AND PROTECTION
TO THE PATIENT.

BENEFITS

- Facilitates the insertion of endoscopic probes.
- Protection against damage to the optical fiber caused by the bite.
- Facilitates oxygen delivery during the procedure.
- Allows the use of the capnograph.

INCLUDED ITEMS

- 01 Orophar - Cannula with double lumen.
- 01 Band for fixation.



INDICATION

- Indicated for airway patency in anesthetic procedures in deep sedation and endoscopic procedures, only in adult patients.

540 GMI DOUBLE LUMEN OROPHARYNGEAL CANNULA - OROPHAR ANVISA Registration No. 80423540061			
Code	SMALLER LUMEN DIMENSION (mm)	LARGER LUMEN DIMENSION (mm)	Cannula Length (mm)
540.5.90	↕ 13x13	Ø 22	90
540.5.100	↕ 13x13	Ø 22	100