

# GMI CATHETER FOR PERITONEAL DIALYSIS

712 - URINARY LINE



## ACCURACY

- Two cuffs that facilitate stabilization.

## BIOCOMPATIBLE

- Developed in medical grade silicone.
- High durability.


## COMFORT

- Extreme proximal rounded.
- Easy insertion.
- Excellent flow.

THE GMI  
 PERITONEAL  
 DIALYSIS  
 CATHETER  
 PROVIDES SAFE  
 AND EASY  
 ACCESS FOR  
 PERFORMING  
 PERITONEAL  
 DIALYSIS.

### PRACTICALITY

- Radiopaque line that facilitates visualization.



### INDICATION

The catheter for peritoneal dialysis with two GMI cuffs is indicated for patients with acute or chronic renal failure undergoing outpatient or home peritoneal dialysis. It provides an easy and safe access to perform peritoneal dialysis, through the peritoneum, as a means for infusion of the dialysis solution.

### CHARACTERISTICS

- 100% medical grade silicone.
- Radiopaque line.
- Sterile.
- Multiperforated with two cuffs.
- PVC extension of 10cm with *luer lock* connection and two adapters.
- Available in neonatal, pediatric and adult sizes.

712   PERITONEAL DIALYSIS CATHETER WITH TWO GMI CUFFS ANVISA No. 80423540036				
Code	FR	EXTERNAL DIAMETER (mm)	INTERNAL DIAMETER (mm)	Length (cm)
712-09-31	09	3.0	1.75	31
712-09-37	09	3.0	1.75	37
712-12-37	12	4.0	2.50	37
712-15-42	15	5.0	3.25	42