


# GMI DOUBLE PIGTAIL URETERAL CATHETER

405/406H - URINARY LINE



## PRACTICALITY

- Easy insertion and adaptation.



## SAFETY

- Material confectioned in polyurethane that guarantees biocompatibility in use.

## VERSATILITY

- - Hydrophilic or teflonized guide that allows a choice suitable for each patient.


THE DOUBLE-J  
CATHETER IS  
USED TO  
PROMOTE FREE  
DRAINAGE OF  
URINE FROM  
THE KIDNEY TO  
THE BLADDER  
IN URINARY  
TRACT  
SURGERY.

### EFFECTIVENESS

- Curved tip, pig tail type, which allows better anchorage of the catheter.
- Multi-hole catheter for excellent drainage.

### INCLUDED ITEMS

- Ureteral Catheter
- Mandrel Guide
- Positioner
- Clamps



### INDICATION

- The Double J Catheter is a flexible tube, which is placed in the ureter during or after surgery to prevent or relieve a blockage. This tube has a 360° bend at both ends, designed to prevent the catheter from moving into the bladder or into the kidney.

GMI DOUBLE PIGTAIL URETERAL CATHETER WITH TEFLONIZED GUIDE						
ANVISA Registration No. 80423540004						
CODE 405A   2 OPEN ENDS						
Guide	FR	LENGTH 20cm	LENGTH 24cm	LENGTH 26cm	LENGTH 28cm	External Diameter
0.032''	4,7	405-47-20A	405-47-24A	405-47-26A	405-47-28A	1,60
0.035''	06	405-06-20A	405-06-24A	405-06-26A	405-06-28A	2,05
0.035''	07	405-07-20A	405-07-24A	405-07-26A	405-07-28A	2,30
0.035''	08	405-08-20A	405-08-24A	405-08-26A	405-08-28A	2,70
CODE 405   1 CLOSED END AND ANOTHER WITH OPEN END						
Guide	FR	LENGTH 20cm	LENGTH 24cm	LENGTH 26cm	LENGTH 28cm	External Diameter
0.032''	4,7	405-47-20	405-47-24	405-47-26	405-47-28	1,60
0.035''	06	405-06-20	405-06-24	405-06-26	405-06-28	2,05
0.035''	07	405-07-20	405-07-24	405-07-26	405-07-28	2,30
0.035''	08	405-08-20	405-08-24	405-08-26	405-08-28	2,70
GMI DOUBLE PIGTAIL URETERAL CATHETER WITH HYDROPHILIC GUIDE						
ANVISA Registration No. °80423540020						
CODE 406H-A   2 OPEN ENDS						
Guide	FR	LENGTH 20cm	LENGTH 24cm	LENGTH 26cm	LENGTH 28cm	External Diameter
0.032''	4,7	406H-47-20A	406H-47-24A	406H-47-26A	406H-47-28A	1,60
0.035''	06	406H-06-20A	406H-06-24A	406H-06-26A	406H-06-28A	2,05
0.035''	07	406H-07-20A	406H-07-24A	406H-07-26A	406H-07-28A	2,30
0.035''	08	406H-08-20A	406H-08-24A	406H-08-26A	406H-08-28A	2,70
CODE 406H   1 CLOSED END AND OPEN END						
Guide	FR	LENGTH 20cm	LENGTH 24cm	LENGTH 26cm	LENGTH 28cm	External Diameter
0.032''	4,7	406H-47-20	406H-47-24	406H-47-26	406H-47-28	1,60
0.035''	06	406H-06-20	406H-06-24	406H-06-26	406H-06-28	2,05
0.035''	07	406H-07-20	406H-07-24	406H-07-26	406H-07-28	2,30
0.035''	08	406H-08-20	406H-08-24	406H-08-26	406H-08-28	2,70