

**BIOCOMPATIBLE**


- 100% medical grade silicone.
- High durability.

**VERSALITY**


- Available in several sizes, in the options 2-way or 3-way

**SAFETY**

- Radiopaque line for easy visualization by image.
- Inflatable balloon for anchoring.



THE  
**GABIFOLLEY, GMI**  
**FOLEY PROBE, IS**  
 AVAILABLE IN  
 SEVERAL SIZES  
 FOR DIFFERENT  
 PATIENT  
 PROFILES AND  
 HAS ALL THE  
 ADVANTAGES OF  
 MEDICAL GRADE  
 SILICONE.



### COMFORT

Malleable and with an atraumatic rounded tip.

### CHARACTERISTICS

Made of 100% sterilized medical grade silicone.

- It has a rounded tip and polished holes for atraumatic passage.
- Nylon guidewire, available in sizes 06, 08 and 10 fr, which make it easier to get through.

### INDICATION

Gabifolley, a GMI silicone foley probe, is indicated for patients who require strict control of urinary output, emptying and bladder irrigation.

450 /475   GMI FOLEY PROBE WITH BALLOON - GABIFOLLEY				
ANVISA Registration No. 80423540007				
TWO-WAY				
Code	FR	External Diameter (mm)	Balloon Capacity (ml)	Length (cm)
450-06	6.0	2.00	1.5	30
450-08	8.0	2.70	3.0	30
450-10	10.0	3.30	3.0	30
450-12	12.0	4.00	10/15*	40
450-14	14.0	4.70	15	40
450-16	16.0	5.30	15	40
450-18	18.0	6.00	15	40
450-20	20.0	6.70	15	40
450-22	22.0	7.30	15	40
450-24	24.0	8.00	30	40
*Available in 10 or 15 ml options, check label.				
THREE-WAY				
Code	FR	External Diameter (mm)	Balloon Capacity (ml)	Length (cm)
475-16	16.0	5.30	30	40
475-18	18.0	6.00	30	40
475-20	20.0	6.70	30	40
475-22	22.0	7.30	30	40
475-24	24.0	8.00	30	40