


BIOCOMPATIBLE

- 100% medical grade silicone.
- High durability.



SAFETY

- Precise and versatile connection.

PRACTICALITY

- Easy and atraumatic insertion.

THE GMI
NELATON
PROBE
PROVIDES
THE
COMPLETE
EMPTYING OF
THE
BLADDER
EFFECTIVELY
AND SAFELY



COMFORT

- Traumatic finish.
- Polished hole at the distal end.

INDICATION

The GMI Nelaton Probe is indicated for immediate bladder emptying and intermittent bladder catheterization. In patients in a hospital or home environment.

CHARACTERISTICS

- 100% medical grade silicone.
- Rounded tip and hole at the distal end.
- Length of 40cm and radiopaque line.
- Available in pediatric and adult neonatal sizes.
- Sterile.

420 GMI NELATON PROBE No. 80423540010			
Code	Internal Diameter (mm)	External Diameter (mm)	FR
420-04	0.75	1.40	04
420-06	1.00	2.00	06
420-08	1.50	2.70	08
420-10	2.00	3.30	10
420-12	2.50	4.00	12
420-14	3.00	4.70	14