


**HYGIENIC**

- Velcro protective flap for fastening reinforcement and protection against dirt.


**FIXING**

- Effective and delicate adhesiveness to the skin.



- Anatomic internal design with a perfect fit for the catheter hub's fins.

GABILOCK  
 PROMOTES  
 PROTECTION  
 AND SAFETY,  
 STABILIZING  
 THE PICC  
 CATHETER,  
 BRINGING  
 COMFORT TO  
 THE PATIENT  
 AND  
 EFFICIENCY IN  
 ITS USE



**PREVENTION**

- Made of hypoallergenic material and adhesive.

**INDICATION**

Gabilock is a fixation adhesive for PICC catheters, which ensures the stabilization of the catheter after puncture. It promotes safety in positioning and during its use, considerably reducing the risk of displacement and loss of the catheter.

AD-650   GABILOCK CATHETER FIXING ADHESIVE			
ANVISA Registration No. 80423540049			
Code	LENGTH (mm)	WIDTH (mm)	Use
AD-650-N	43,18 +/- 1,0	15,24 +/- 1,0	Neonatal
AD650-PA	56,92 +/- 1,0	15,24 +/- 1,0	Pediatric Adult