


GMI DOUBLE LUMEN PICC CATHETER

581 - VASCULAR LINE



DIVERSITY

- Two lumens that allow multiple infusions.

COMFORT

- Anatomical shape of the hub providing safe stabilization.

PRACTICALITY

- Made of polyurethane, material that facilitates the passage of the catheter.

THE PICC
DOUBLE
LUMEN
CATHETER IS
INSERTED
THROUGH A
PERIPHERAL
PUNCTURE
INTO A
CENTRAL
ACCESS IN
PATIENTS WHO
REQUIRE
MEDIUM- AND
LONG-TERM
INFUSION.

SAFETY

- Centimetered catheter that facilitates positioning during passage.
- Impregnated in barium sulfate, facilitating its visualization.



INDICATION

Catheter of choice for neonatal use, whenever there is a need for long-term intravenous therapy.

It has several advantages such as greater comfort and safety for the patient, low risk of complications and associated infections, and cost reduction for the institution. Insertion by inserter into peripheral vein reaching the central position.

It allows the infusion of several intravenous therapies such as TPN, chemotherapy, hyperosmolar, vesicant or irritant, with the safety of simultaneous use of more than one drug, avoiding drug interaction.

INCLUDED ITEMS

- 01 Catheter 100% Polyurethane
- 01 Peel Off Introducer
- 01 Tape measure
- 01 Tweezers
- 01 Tourniquet
- 01 Plastic Sheath
- 01 Attachment Sticker
- 01 Transparent Film Coverage

581 GMI POLYURETHANE DOUBLE LUMEN PICC CATHETER					
ANVISA Registration No. 80423540078					
Code	FR	NUMBER OF WAYS	TOTAL PRIME (ml)	MAXIMUM PRESSURE (PSI)	Length (cm)
581-02-30-2	2,0	02	0,32	220	30